

Oxy Foam




Foaming Cup



Power Foamer

**Exclusive
Foaming Systems**

Quality of a wet cleaning foam



The % moisture of a foam is directly linked to its cleaning power and its capacity to remove soiling in all the micro-recesses of the skin.

Oxy Foam



Active foam for teat hygiene

Oxy Foam




Well proven foaming technology, easy to use,
safe and economical


An active and wet foam
ready in seconds with the
exclusive Ecolab foaming
cup or Power Foamer.




Fill the teat cup refill
bottle to 2/3 full



Shake to and fro



Squeeze the bottle and push the foam
up until overflows the teat cup



Or squeeze simply on the trigger

Easy to apply



Apply the foam to the teats



Allow sufficient contact time



Wipe off with one-use paper



The results : clean teats

Oxy Foam, an active foam with following benefits :

- The wet foam completely covers the teats for an optimal cleaning.
 - **Oxy Foam** becomes directly active on contact and blocks cross contamination.
 - **Oxy Foam** is hypo allergenic and softens and conditions skin thanks to its cosmetic ingredients.
 - **Oxy Foam** is highly visible on the teats and is easily wiped off before milking.
 - More effective application using less product. No product loss from "run off".
- Oxy Foam** is a formulation that meets the needs of modern food safety requirements.

The foaming components in **Oxy Foam** are specially designed to generate a wet foam with enhanced technical properties when used either with the exclusive Ecolab foaming cup or Power Foamer. The constant micro-bubbling of the foam together with its adhesive consistency provides complete removal of soiling from the teats as well as treating sensitive areas of the skin.